




# **Standard electrode positions and reference systems**

Nico Franzmeier




# EEG - overview

- What is recorded?
  - Electrical activity of the cortex
  - Voltage fluctuations from ionic current within neurons




# EEG - challenges

- In order to be **comparable** between subjects and sites and to fulfill clinical and scientific requirements, the EEG signal should be **recorded in a standardized manner**
- **Science:** Different brain regions are involved in different cognitive processes



# EEG - challenges

- In order to be **comparable** between subjects and sites and to fulfill clinical or scientific requirements, the EEG signal should be **recorded in a standardized manner**
- **Clinic:** Different brain regions are involved in different etiologies of neurological disorders (epilepsy → epilepsy surgery)




# EEG - challenges

- A standardized EEG electrode positioning system is required




# The 10-20 System of Head Measurement and Electrode Placement

- Standardized placement system of EEG electrodes established in 1958 by H. Jasper, a canadian psychologist
- Meant to facilitate communication and comparison of results
- Provides for adequate coverage of all parts of the head/brain



# The 10-20 System of Head Measurement and Electrode Placement


- The numbers '10' and '20' refer to the distances between adjacent electrodes
- Usually, 21 (19 active + ground & reference) electrodes are used, but others can be added to improve spatial resolution



# The 10-20 System of Head Measurement and Electrode Placement

## Terminology:



- Each electrode has a **letter** to identify the lobe and a **number** to identify the hemisphere location
- **Even numbers** (2,4,6,7) refer to electrodes on the right hemisphere
- **Odd numbers** (1,3,5,7,) refer to electrodes on the left hemisphere
- Electrodes along the mid sagittal plane have the **subscript “z”**



Electrode	Lobe
F	Frontal
T	Temporal
C	Central *
P	Parietal
O	Occipital


# The 10-20 System: Basic steps

- Step 1: Four anatomical landmarks are used to position the electrodes





# The 10-20 System: Basic steps

- Step 1: Measure over the center line of the scalp, from the Nasion (bridge of the nose) to the Inion (external occipital protuberance)





# The 10-20 System: Basic steps

- Step 2: Measure and mark 50% of your total. This is the preliminary Cz mark




# The 10-20 System: Basic steps

- Step 3: Measure and mark 10% up from the Nasion and 10% up from the Inion. These are the preliminary Fpz and Oz marks



# The 10-20 System: Basic steps


- Today EEG-caps are used commonly, thus placing individual electrodes is done only in exceptional cases




=> What if the spatial resolution is not high enough for my purposes?

# Extending the 10-20 system


- In a standard 10-20 setup, 21 electrodes are used
- Other electrodes can be added, to
  - Increase spatial resolution (> 256 electrodes are possible)
  - Record from specific areas of interest
  - Monitor other electrical activity (e.g. ECG, EMG, eye movement)



10-5



10-10



10-20

# Extending the 20-10 system

- How does a higher number of electrodes improve spatial resolution?
- Mean interelectrode distance  
(Gevins, 1994)
  - Classical 10-20 system: 6 cm
  - ...64 electrodes: 3.3 cm
  - ...128 electrodes: 2.25 cm
  - ...256 electrodes: 1.6 cm





124-channel




256-channel

# If more than 256 electrodes still aren't enough



# Standardized electrode positioning and underlying brain regions



- ⇒ Clinical purpose: Detect the origin of a seizure
- ⇒ Scientific purpose: Investigate brain regions/networks and their involvement in cognitive processes/neurological disorders

# The EEG signal

- Electric potentials are only defined with respect to a reference electrode, i.e. an arbitrarily chosen “zero level”
- The EEG signal represents a difference in voltage between an active electrode and a reference
- The choice of the reference may depend on the purpose of the recording

# EEG Referencing

- Ideally, the reference
  - would be affected by global voltage changes in the same manner as all other electrodes
    - Unspecific activity is subtracted out by referencing (e.g. slow voltage shifts due to sweating)
  - should not pick up signals that are not intended to be recorded (i.e. ECG)
    - Hence It should be placed in an presumed “inactive” zone

# Referencing montages in EEG

## ■ Common reference


- All electrodes use the same reference (e.g. mastoid, earlobe, Cz, Nose)
- **But:** unilateral Reference electrodes can introduce asymmetries

## ■ Average reference

- Average between two electrodes (e.g. left and right earlobe or mastoids)
- Common average reference: Global mean of all electrodes
- Local average reference: Mean of the 4 or 8 surrounding electrodes

# Common referencing methods in EEG

- In Addition to the active and reference electrode, typical EEG systems also use a ground electrode



Voltage is measured between ACTIVE and GROUND (A - G)

Voltage is measured between REFERENCE and GROUND (R - G)

Output is difference between these voltages

$$(A - G) - (R - G) = A - R$$

It's as if the ground does not exist

Any noise in common to A and R will be eliminated



# Thanks for your attention

